

Want More Options To Restore PREVAIL® 4/3 Implants? Try New BIOMET 3i 3.4mm Restorative Components Today!



INTRODUCING NEW BIOMET 3i 3.4mm Restorative Components



**Straight And Pre-Angled
GingiHue® Abutments For
Stock Restorative Options**

**Immediately Provisionalize
With Straight And Pre-Angled
PreFormance® Posts**



**Soft-Tissue Healing With
New Healing Abutments**

**Twist Lock™ Impression
Copings Now Available!**



BIOMET 3i 3.4mm Restorative Components Provide **Clinicians** One Solution At A Time

- Restore PREVAIL® 4/3 Implants With New 3.4mm GingiHue® Abutments
- New 3.4mm PreFormance® Posts For Immediate Aesthetic Options
- 3.4mm Healing Abutments Emerge To 5.0mm For Soft-Tissue Healing
- Impression Copings Emerge To A 5.0mm Emergence Profile (EP®) For Closed Tray Impression Techniques

BIOMET 3i
PROVIDING SOLUTIONS - ONE PATIENT AT A TIME™

Restore PREVAIL® 4/3 Implants Today With New 3.4mm Restorative Components!

Item	Description
GingiHue® Abutments	
IMAP352G	Certain® GingiHue® Abutment 3.4MM(D) X 5.0MM(P) X 2.0MM(H)
IMAP354G	Certain® GingiHue® Abutment 3.4MM(D) X 5.0MM(P) X 4.0MM(H)
IMPAP352G	Certain® GingiHue® 15° Pre-Angled Abutment 3.4MM(D) X 5.0MM(P) X 2.0MM(H)
IMPAP354G	Certain® GingiHue® 15° Pre-Angled Abutment 3.4MM(D) X 5.0MM(P) X 4.0MM(H)
Emergence Profile (EP®) One-Piece Healing Abutments	
IMHA352	Certain® EP® Healing Abutment 3.4MM(D) X 5.0MM(P) X 2.0MM(H)
IMHA353	Certain® EP® Healing Abutment 3.4MM(D) X 5.0MM(P) X 3.0MM(H)
IMHA354	Certain® EP® Healing Abutment 3.4MM(D) X 5.0MM(P) X 4.0MM(H)
IMHA356	Certain® EP® Healing Abutment 3.4MM(D) X 5.0MM(P) X 6.0MM(H)
Pick-Up Impression & Twist Lock™ Impression Copings	
IMIC35	Certain® EP® Pick-Up Coping 3.4MM(D) X 5.0MM(P)
IMIT33	Certain® EP® Twist Lock™ Coping 3.4MM(D) X 3.8MM(P)
IMIT35	Certain® EP® Twist Lock™ Coping 3.4MM(D) X 5.0MM(P)
PreFormance® Posts	
IMPAPF354	Certain® PreFormance® 15° Pre-Angled Post 3.4MM(D) X 5.0MM(P) X 4.0MM(H)
IMPAPF356	Certain® PreFormance® 15° Pre-Angled Post 3.4MM(D) X 5.0MM(P) X 6.0MM(H)
IMPFP354	Certain® PreFormance® Post 3.4MM(D) X 5.0MM(P) X 4.0MM(H)
IMPFP356	Certain® PreFormance® Post 3.4MM(D) X 5.0MM(P) X 6.0MM(H)
Manual Bone Profiler	
IMMBP35	Certain® Manual Bone Profiler 3.4MM(D) X 5.0MM(P)

(D = Diameter, P = Platform, H = Height)

BIOMET 3i
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Outside The U.S.: +1-561-776-6700
Fax: +1-561-776-1272
www.biomet3i.com

EC REP **BIOMET 3i**
Europe, Middle East & Africa
WTC Almeda Park, Ed. 1, Planta 1ª
Pl. de la Pau, s/n
08940, Cornellà de Llobregat
(Barcelona) Spain
Phone: +34-93-470-55-00
Fax: +34-93-371-78-49

Enhanced
Knowledge



Comprehensive
Treatment

Business
Excellence



ART1174

REV A 11/11

Not Available In All Markets. Please Consult Your Local BIOMET 3i Sales Representative For Availability.

BIOMET 3i™
PROVIDING SOLUTIONS – ONE PATIENT AT A TIME™



Join
Us



Follow
Us



Watch
Us



Download
It